






AlfaGate
IMPLANTS YOU CAN TRUST

Alfa-Lock Attachment System






Low Profile Attachment System
The Best Choice For Implant Retained over-dentures.

- 
- 
- 
- 
- 01 SP Platform
 - 02 S-Line Platform
 - 03 NP Platform
 - 04 RP Platform
 - 05 RC Platform
 - 06 Replacment&Tools



Alfa-Lock Attachment



- Additional Info: Recommended torque **25Ncm**. Using hex driver 1.25mm.
- Material: Titanium

Height mm	Ref.
1mm	AGM-106-1
2mm	AGM-106-2
3mm	AGM-106-3
4mm	AGM-106-4
5mm	AGM-106-5

Alfa-Lock Set For SP Connection



Height mm	Ref.
1mm	AGM-106-1K
2mm	AGM-106-2K
3mm	AGM-106-3K
4mm	AGM-106-4K
5mm	AGM-106-5K


All Alfa-Lock Sets include:

Description	Hard plastic cap	Standard plastic cap	Soft plastic cap
Ref.	AGM-106-PH 	AGM-106-P 	AGM-106-PS 
Metal housing		Protective disk	
AGM-106-M		AGM-106-D 	

*Alfa-Lock can be sold separately.



Alfa-Lock Attachment



- Additional Info: Recommended torque **20 Ncm**. Using hex driver 1.25 mm.
- Material: Titanium

Height mm	Ref.
1mm	AGM-106-1S
2mm	AGM-106-2S
3mm	AGM-106-3S
4mm	AGM-106-4S
5mm	AGM-106-5S

Alfa-Lock Set For S-Line Connection



Height mm	Ref.
1mm	AGM-106-1SK
2mm	AGM-106-2SK
3mm	AGM-106-3SK
4mm	AGM-106-4SK
5mm	AGM-106-5SK


All Alfa-Lock Sets include:

Description Ref.	Hard plastic cap AGM-106-PH	Standard plastic cap AGM-106-P	Soft plastic cap AGM-106-PS
Metal housing AGM-106-M	Protective disk AGM-106-D		

*Alfa-Lock can be sold separately.



Alfa-Lock Attachment



- Additional Info: Recommended torque **20 Ncm**. Using hex driver 1.25 mm.
- Material: Titanium

Height mm	Ref.
1mm	AGM-106-1NP
2mm	AGM-106-2NP
3mm	AGM-106-3NP
4mm	AGM-106-4NP
5mm	AGM-106-5NP

Alfa-Lock Set For NP Connection



Height mm	Ref.
1mm	AGM-106-1NPK
2mm	AGM-106-2NPK
3mm	AGM-106-3NPK
4mm	AGM-106-4NPK
5mm	AGM-106-5NPK


All Alfa-Lock Sets include:

Description	Hard plastic cap		Standard plastic cap		Soft plastic cap	
Ref.	AGM-106-PH		AGM-106-P		AGM-106-PS	
Metal housing		Protective disk				
AGM-106-M		AGM-106-D				

*Alfa-Lock can be sold separately.



Alfa-Lock Attachment



- Additional Info: Recommended torque **25Ncm**. Using hex driver 1.25mm.
- Material: Titanium

Height mm	Ref.
1mm	AGM-106-1RP
2mm	AGM-106-2RP
3mm	AGM-106-3RP
4mm	AGM-106-4RP
5mm	AGM-106-5RP

Alfa-Lock Set For RP Connection



Height mm	Ref.
1mm	AGM-106-1RPK
2mm	AGM-106-2RPK
3mm	AGM-106-3RPK
4mm	AGM-106-4RPK
5mm	AGM-106-5RPK


All Alfa-Lock Sets include:

Description Ref.	Hard plastic cap AGM-106-PH	Standard plastic cap AGM-106-P	Soft plastic cap AGM-106-PS
Metal housing AGM-106-M	Protective disk AGM-106-D		

*Alfa-Lock can be sold separately.



Alfa-Lock Attachment



- Additional Info: Recommended torque **25Ncm**. Using hex driver 1.2mm.
- Material: Titanium

Height mm	Ref.
1mm	AGM-106-1RC
2mm	AGM-106-2RC
3mm	AGM-106-3RC
4mm	AGM-106-4RC
5mm	AGM-106-5RC

Alfa-Lock Set For RCConnection



Height mm	Ref.
1mm	AGM-106-1RCK
2mm	AGM-106-2RCK
3mm	AGM-106-3RCK
4mm	AGM-106-4RCK
5mm	AGM-106-5RCK

All Alfa-Lock Sets include:

Description Ref.	Hard plastic cap AGM-106-PH		Standard plastic cap AGM-106-P		Soft plastic cap AGM-106-PS	
Metal housing AGM-106-M		Protective disk AGM-106-D				

***Alfa-Lock can be sold separately.**

Alfa-Lock Replacement Set

Ref. AGM-106-1T


Description	Standard plastic cap		Soft plastic cap	
Ref.	AGM-106-P		AGM-106-PS	
Hard plastic cap			Metal housing	
AGM-106-PH		AGM-106-M		Protective disk
			AGM-106-D	




Alfa-Lock Multi Tool


Ref. IT-AL

For inserion /Extraction Plastic caps.





Alfa Gate has 25 years of global success in manufacturing and supplying high-quality dental implants and prosthetic solutions curated to all dentists. The effectiveness and superior performance of our products are a result of successful R&D that ensures the greatest quality.



w w w . a l f a - g a t e . c o m

